015-21001

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## 1.1 preface

**Please read this instruction manual carefully before using this product.**

Thank you for purchasing the GC-400 portable gas chromatograph VOC special model (this product) from Shanghai Huixian Instruments Co., Ltd.

This manual describes the operation method of this product (hardware installation and software analysis application). Please read this instruction carefully and use the product correctly before using this product.

Please keep this instruction manual properly for future reference.

The content recorded in this manual is based that the user has the basic operation knowledge of Windows.

This manual is a double edition in Chinese and English

## material particulars

- If the user or place changes, please transfer this instruction to subsequent users.
- If this manual is lost or damaged, please contact your dealer in each area or at our after-sales department immediately.
- To ensure safe operation, please contact the professional after-sales engineer of the company for installation, commissioning or maintenance.

## statement

This instruction manual is subject to changes without notice.

This manual strives to be accurate, please have any mistakes or omissions, please forgive.

The copyright of this manual belongs to Shanghai Huixian Instrument and Equipment Co., Ltd. Part or all shall not be reproduced without the permission of the Company.

Microsoft, Windows and Windows 10 are registered trademarks of American Microsoft Corporation in the United States and other countries.

The Adobe, Adobe logo and Adobe Reader are registered trademarks of American Adobe Systems Incorporated in the United States and other countries.

Other company names and product names recorded in this instruction are the trademarks and registered trademarks of each Company. In addition, the TM, markers are not clearly specified in this specification.

**The Microsoft Windows 10 operating system is called the "Windows 10".**

The supply period of spare parts of this product is seven years after the production stop. Please understand that spare parts may not be available later. However, the supply period of spare parts please refer to the provisions of relevant manufacturers

## 1.2 guarantee

The warranty content of this product is as follows:

### 1. defects liability period

No special extended warranty contract for this product is provided within 12 months of the date of equipment installation or within 13 months after arrival.

### 2. warranty content

Failure caused by the Company will be repaired or replaced with parts (including software) during the warranty. However, the same model may not be available for short-lived products such as personal computers, their peripherals and components.

### 3. scope cover

- (1) In no event shall the Company be liable to any user for missed work costs, indirect damages and derivative damages. Nor is it liable for damages brought to the third user.
- (2) The maximum amount of the Company shall be limited to the factory price or sales price.

### 4. waiver

The following faults are not under the warranty:

- 1) misoperation.

Repair or modification of this product by 2) not Company or other companies designated by the Company.

- 3) is used with non-hardware or software designated by the company.

4) This product failure due to computer viruses and software and data damage including basic software.

5) faults caused by power failure such as power failure or sudden voltage reduction and software and data damage including basic software.

Failure caused by 6) error shutdown and software and data damage, including basic software.

- 7) is not due to the product itself.

8) fault caused by the use of this product in harsh environments such as high temperature and humidity, corrosive gas or vibration.

9) has faults caused by fires, earthquakes, other natural disasters, contamination of radioactive and harmful substances, as well as irresistible accidents such as war, riots and crime.

Failure caused by 10) when dismantling or transporting products after installation.

11) consumables or equivalent components.

### 1.3 Warning ID

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Common labels in this

**warn**

For situations that can cause death or serious

**take**

For situations that may cause mild injury or

**pay**

Provide additional information to ensure

## 2.1 Product Introduction:

**Type GC-400 portable VOC analyzer adopts the separation of methane + FID detection technology by double quantitative cycle quantification and column separation of methane, total hydrocarbon and non-methane total hydrocarbon concentration.**

### **Main characteristics of type GC-400 portable VOC total hydrocarbon analyzer:**

Using the chromatographic column separation method + FID detection technique, the high separation efficiency, non-methane total hydrocarbon concentration analysis period was 60 seconds;

Automatic determination of methane and non-methane total hydrocarbon: fully automatic flow path design, single valve and double channel sample inlet, double quantitative ring, and automatic determination of methane and non-methane total hydrocarbon values at one time;

V O C special station of win10 operating system, automatic control of measurement process; automatic conversion of ppm mg/m<sup>3</sup> unit without manual calculation.

Support for bulk data processing reports;

Built-in printers support automatic summary of data from field printing.

The carrier uses nitrogen with little oxygen interference; and the software has automatic oxygen deduction function;

Using low-pressure solid-state hydrogen storage bottle, safe and reliable, easy to carry and transport, can be used according to the set flow of more than 9 hours;

The whole high temperature heat companion, effectively solve the problem of condensation loss of samples under high temperature and high humidity gas occasions;

The supporting heat companion pipe sampling gun is a polytetrafluoride pipeline, with little adsorption impact on VOC volatilization interference detection;

Built-in lithium battery, can be preheated up in advance, save waiting time, to facilitate the use of equipment transfer;

Built-in large-capacity lithium battery and thermal companion gun can work for 4 hours, and choose with a large battery life of up to 15 hours;

Optional working condition gun, Bluetooth connection measurement and calculation includes dynamic pressure, static pressure, full pressure, flue gas flow rate, flue gas temperature, moisture content, conversion concentration, oxygen content and other parameters;

Built-in printer, detection data field printing, and support U disk export;

The whole machine is less than 15kg, portable and small, convenient for on-site use.


### **Main technical parameters:**

|                          |                                                                                                                        |
|--------------------------|------------------------------------------------------------------------------------------------------------------------|
| function                 | Analysis of methane; total hydrocarbon; non-methane total hydrocarbon; benzene series and other V O C s (are optional) |
| size                     | 514 mm ( L ) x 340 mm ( W ) x 150 mm ( H )                                                                             |
| Host weight              | ≤13 K g                                                                                                                |
| monitor                  | FID detector; meet 10 <sup>6</sup> Dynamic linear range                                                                |
| Principle of detection   | GC-FID; chromatcolumn separation                                                                                       |
| Traffic control accuracy | 0.001kpa                                                                                                               |
| Analysis cycle           | NMHC:60s                                                                                                               |
| range ability            | 0 ~ 50,000 p p m; Adjustable                                                                                           |
| Minimum detection limit  | ≤ 0.03mg/m <sup>3</sup>                                                                                                |

|                       |                                                                      |
|-----------------------|----------------------------------------------------------------------|
| repetitiveness        | $RSD \leq 0.5\%$                                                     |
| accuracy              | $\leq \text{range} \pm 1\%$                                          |
| operation panel       | <b>touch control</b>                                                 |
| communication         | GPS positioning (optional), printing, data upload, WIFI, Bluetooth   |
| source                | Built-in power supply, with a battery life of more than 4 hours      |
| Probe temperature     | The constant temperature of $\sim 200^{\circ}\text{C}$ is adjustable |
| Valve box temperature | $200^{\circ}\text{C}$ Max .                                          |
| ambient temperature   | $-15^{\circ}\text{C} \sim 50^{\circ}\text{C}$                        |
| work environment      | 5 - 90 % R H                                                         |


## 2.2 operational principle

The GC-400 portable VOC instrument is composed of gas circuit system, sample separation system, detector system, signal processing system, data processing system, power supply system, temperature control unit, display circuit, and analysis software. The samples are separated from the gas belt into the chromatography column and then converted to electrical signal by the detector. It is amplified and plastic by the data acquisition circuit, transforms analog amount to digital amount, and the internal calibration data is transferred for calculation and processing, and then the measured value is displayed and printed by the display circuit.

Analyze the flow chart

## 2.3 Introduction of functional area of equipment panel

graph 1

| number | explain                                                                             |
|--------|-------------------------------------------------------------------------------------|
| 1      | Pressure adjustment area for carrying gas and hydrogen air pressure flow adjustment |
| 2      | Pressure display area                                                               |
| 3      | Auxiliary heating area, used for solid-state hydrogen storage gas heating           |
| 4      | Gas link mouth                                                                      |
| 5      | Sampling gun interface for sampling gun control and gas collection                  |
| 6      | printer                                                                             |
| 7      | Battery, used for battery replacement and equipment maintenance                     |
| 8      | mains switch                                                                        |
| 9      | External power supply and charging interface                                        |
| 10     | radiator fan                                                                        |
| 11     | USB joggle                                                                          |
| 12     | display control unit                                                                |

### **3.1 Use environment requirements**

Consider the following points to ensure the safe and correct operation of the equipment.

#### **3.1.1 ambient temperature and humidity**

Operating this product for optimal performance at the temperature and humidity in the following range

temperature range: 18℃-28℃

Relative humidity range: 50%-60%

operating temperature range: 5℃-40℃

Operating humidity range: 5%-90%

#### **3.1.2 placement platform requirements**

Place the instrument in firm, stable, and flat places

(GC-400 Portable Gas chromatograph VOC weighs about 13kg)

#### **3.1.3 corrosive gas and dust**

Avoid exposure to corrosive gases and excessive dust to extend the service life of the equipment and maintain its optimal performance.

#### **3.1.4 电 magnetic field and power supply device noise**

This instrument cannot be used near a strong magnetic field. The noise generated by the power supply device must be very low or minimal. These factors have the potential to cause an instrument failure.

#### **3.1.5 other considerations for 3.1.5**

Avoid the following conditions to maintain optimal performance:

3.1.5.1) ambient temperature fluctuations.

3.1.5.2) heating or air conditioning changes the temperature.

The 3.1.5.3) is directly exposed to sunlight.

The 3.1.5.4) will vibrate strongly

### 3.2 shift

Be careful to move this product without collision or vibration.

3.2.1. weighs approximately 13 kg (type GC-400 V).

3.2.2. must be lightly.

3.2.3. shall not be inverted when handling.

## 4 install

### 4.1.1 electrical device and lines

**Warn**

**Pay attention to the high-voltage power when using the external power supply!**

A Turn off the power supply on the wiring board before connecting the wire to the wiring board.

The b power supply device must have a circuit breaker

C Do not place heavy objects on the wire

### 4.1.2 electric device voltage

Use the power supply as described below to maintain the optimal performance of the instrument.

Recommended voltage using power supply: lithium battery charger 25.2V 5A

Frequency of 50/ 60Hz

Built-in polymer lithium battery charging voltage: 25.2V 5A

Frequency of 50/ 60Hz

### 4.1.3 electric capacity

Model number: 24V30000

Nominal voltage: DC24V

Operating voltage: 19.2-25.2V

capacity: 30000MA/H

Cell: lithium polymer battery


Size: 147mm\*110mm\*62mm

weight: 750g

Sustainable maximum current: 10A

Overcurrent protection value: 50A

pay attention to Voltage fluctuation of power supply device fluctuates or insufficient capacity will affect instrument performance.

## 4.2 Gas supply and installation

**take care**

4. No gas supply pressure above the maximum pressure shown in the table below will damage the pressure control valve or other components.
5. When sharing with other instruments-one gas supply device, check the operating instructions for all devices including instruments and gas supply,
6. All devices operate within the required pressure range.

The gas and purity values required below can maintain the optimal performance of the instrument.  
FID and TCD except the FID-detector.

### 4.2.1 body type

Type of adenyated gas

(FI filled column D analysis)

Helium or nitrogen can be used. Nitrogen is cheaper in price. (Filled column TCD analysis)

Other substances can be analyzed as carriers using helium or hydrogen at high sensitivity. Because hydrogen is flammable substances, out of safety

Sexual consideration, -like the use of helium gas. However, for the analysis of helium and hydrogen, nitrogen and argon is required as a carrier. It would be convenient to analyze trace fractions in the air using nitrogen, because its nitrogen peak is not detected by the detector when loaded as the gas.

(Capillary column analysis)

Helium gas is more well suited for separation.

Although nitrogen is more ideal than helium, optimal separation conditions are not easily achieved.

Type of caudal blow

(Capillary Column FID)

Both helium and nitrogen can be used. Nitrogen sensitivity is high. (Capillary Column TCD)

The same gas was used for tail blowing and reference gas for capillary TCD analysis. Select the gas type with the same method for carrying the gas.

2. gas purity

Helium gas (gas carrier, tail blowing) : 99.995% or higher

Nitrogen (gas load, tail blowing) : 99.995% or higher

Argon (gas load, tail blowing) : 99.995% or higher

Hydrogen gas (FID detector gas) : 99.995% or higher

Air (FID detector gas) : Dry air (remove oil and other organic ingredients)

Compressed air (which must be compressed and dehumidified with an oil-free compressor)

### 4.2.2 body supply pressure

carrier gas 300-980kPa (hydrogen: 300-500kPa)

Tail blow 300-980kPa

hydrogen 300 - 500kPa

air 400 - 500kPa

Note that the conversion relationship between kPa and Ba is described below


100kPa = 1 bar

The conversion relationship between kPa and kgf/cm<sup>2</sup> is as follows

1 kPa = 1.02 × 10<sup>-2</sup> kgf/cm<sup>2</sup>

The conversion relation of 1 kgf/cm<sup>2</sup> = 98.1 kPa kPa and psi is as follows  
1 kPa = 1.45 × 10<sup>-1</sup> psi  
1 psi = 6.89 kPa

#### 4.3 Notes for high-pressure gas tank

Common considerations are described below.


Reference to national and local special adjustments.

Place the tank vertically in a flat position, not upside down.

The air tank temperature can not exceed 40 ° C. Flammable substances must be placed two meters outside the air tank.

When using high pressure gas and daily leakage inspection, good ventilation must be strictly maintained, especially when using flammable gas (such as hydrogen), equipment must not smoke or ignite open fire must be fire extinguisher equipped with fire extinguisher tools to fix the gas tank to avoid tipping or overturning the pressure valve can only use oil, the main valve of the gas tank must be tightened immediately.

#### 4.3.1 hydrogen

#### 4.4 Connect of gas road supply pipeline


There are two types of link tubes (type P) PTFE and (type M) metal stainless steel or purple copper in this product

##### 4.4.1 气 Road supply pipe specification

1 PTFE  $\Phi 3\text{mm} \times 2\text{mm} \times 1000\text{mm}$

2 Stainless steel  $\Phi 3\text{mm} \times 2\text{mm} \times 1000\text{mm}$

##### Installation of the 4.4.2 减 pressure valve

1 High pressure cylinder 2 pressure gauge 4 rotary connector

If the cylinder-type gas source is used, the pressure reduction valve installation steps are as follows:

- 4.4.2.1 shall respectively unscrew the low pressure outlet heads of two oxygen pressure reduction valves and one hydrogen pressure reduction valve (such as using solid state hydrogen without installing the pressure reduction valve), connect to the pressure reduction valve connector, and turn on the low pressure output control rod (do not tighten);
- 4.4.2.2. Install the pressure relief valve on the cylinder, tighten the screw cap, open the cylinder high pressure valve, the high pressure gauge of the pressure relief valve shall be indicated;
- 4.4.2.3. After closing the cylinder high pressure valve, the high pressure gauge of the pressure reduction valve should not drop, otherwise there is air leakage, should be discharged Only can be use.

## Pipe line between

### 4.4.3 tank and portable PLC

Gas chromatograph connecting pipe

The upper panel of the instrument provides connections to the external pipe. The label are shown below.

Air load, input CARRIER ("L" and "R" refer to the left or right of the two-stream feed.)

make up gas MAKE UP if available

Hydrogen, input HYDROGEN

Air, input AIR


Fig. 1-1 Schematic diagram of the connection between gas cylinder and equipment

GC 400 Portable gas chromatograph external gas pipeline, using polyethylene pipe with outer diameter 30.5 (wall thickness), 20.5 stainless steel liner, O ring, silicone gasket, M81 threaded sealing nut, M81 threaded joint for connection.

The connection method is shown in Figures 2-1.


Fig. 2-1 Schematic diagram of polyethylene pipe gas road connection

The external gas pipeline for the GC 400 portable gas chromatograph can also be connected with a 30.5 outer diameter stainless steel pipe or purple copper pipe, O-ring, silicone gasket, M81 threaded sealing nut, M81 threaded connector.

The connection method is shown in Figures 2-2.


Fig. 2-2 Schematic diagram of metal pipe gas pipeline connection

The connection step are described as follows (take connecting gas source cylinder as an example):

8. Cut off a polyethylene tube of the appropriate length and insert a 20.5 stainless steel liner at each end.
9. Install the M81 sealing nut, phosphocopper rings and 2 O rings into one end of the polyethylene tube.
10. Spin the M81 sealing nut on the cylinder connector (M81) and tighten to ensure a good seal.
11. Install the M81 sealing nut, phosphocopper rings and 2 O rings into the other end of the polyethylene tube.
12. Spin the M81 sealing nut on the corresponding connector (M81) of the purifier, and tighten it to ensure a good seal.
13. The connection steps for other gas pipelines are the same as those described above.
14. If metal pipes such as stainless steel or purple copper pipes are used, the 20.5 stainless steel liner need not be used, and tighten the same as using silicone pad as step 2.3.4.5.

▲ pay attention to:

Enter the pressure requirements to the GC 400 portable gas chromatograph:

4. The gas carrying pressure must be in the 0.343Mpa — 0.392Mpa range (equivalent to 3.5kg / cm<sup>2</sup>—4kg/cm<sup>2</sup>) 。
5. The air pressure must be in the 0.4Mpa — 0.5Mpa range (equivalent to 3kg /


cm<sup>2</sup>—4kg/cm<sup>2</sup>).

6. The pressure of hydrogen must be in the 0.296Mpa —0.392Mpa range (equivalent to 2kg / cm<sup>2</sup>—4kg/cm<sup>2</sup>).

## 5 Inflation and use methods of small gas bottle tanks:

### 5.1 Nitrogen inflation method

Use the QF-2 gas cylinder bridge connector (HS2155-QF0271000) to first connect one end of the inflatable cylinder and tighten, slowly open the main valve slowly to flow out of air from the connector, then connect the inflatable tank and tighten and fully open the main valve of the inflatable cylinder, first close the main valve for 2-5 seconds. Remove the bridge connector to air.


Link Schematic

warn!The maximum pressure resistance of the inflatable tank is 15mpa, and the internal pressure of the inflatable cylinder should not be higher than 15mpa to avoid danger.

It was tested every 5 years

### 5.2 Air inflation method

Use the QF-2 gas cylinder bridge connector (HS2155-QF0271000) to first connect one end of the inflatable cylinder and tighten, slowly open the main valve slowly to flow out of air from the connector, then connect the inflatable tank and tighten and fully open the main valve of the inflatable cylinder, first close the main valve for 2-5 seconds. Remove the bridge connector to air.


#### Link Schematic

warn!The maximum pressure resistance of the inflatable tank is 15mpa, and the internal pressure of the inflatable cylinder should not be higher than 15mpa to avoid danger.

It was tested every 5 years

### 5.3 Description of metal hydrogen storage and use of charging and discharge method

Schematic diagram of the shape and link of hydrogen storage of 5.3.1 metal hydride

#### 5.3.2 充 gas operation process

The hydrogen reservoir is generally filled with a certain amount of argon inside the factory. When first use, argon needs to be pumped out by vacuum pump, or by multiple hydrogen discharge method. A 3mm stainless steel pipe (HS2115-0030021000) is recommended to connect hydrogen storage valves and gas equipment. In order to ensure the purity of hydrogen use, the hydrogen storage device end should be tightened first, open the valve to wash the pipeline, and then tighten the hydrogen equipment end.

**The hydrogen filling process of the 5.3.3 本 instrument is as follows:**

5.3.3.1. hydrogen source: Recommended hydrogen purity of the hydrogen source: > 99.99%.

5.3.3.2. connection: Connect the hydrogen reducer to the hydrogen bottle; connect the hydrogen reducer to the metal hydrogen reservoir with 3 stainless steel pipe. The pipe must be cleaned with hydrogen before the connection to prevent air from the pipe from entering the system. The pipeline cleaning method is: loosen the inlet and outlet trap of the metal hydrogen storage device, adjust the pressure of the hydrogen reducer to 0.1~0.3MPa, let the airflow through the pipeline into the air, to clean the pipeline, row continuously for 5~10 seconds, and then tighten the sleeve and seal the system.

5.3.3.3. leakage detection: Check the interface with soapy water to ensure no air leakage.

5.3.3.4. hydrogen charge: adjust the hydrogen pressure reducer to make the hydrogen pressure reach the required value: open the metal hydride hydrogen storage valve, and start the hydrogen charge. Since the temperature of the metal hydride hydrogen reservoir will increase during

inflation, it is recommended to put the metal hydride hydrogen reservoir in a water bath to cool it.

Recommended hydrogen filling conditions: pressure 0.4~1.0MPa

Time 50~220min (adjusted to different models)

temperature 25℃

#### **5.3.4, hydrogenation operation process**

5.3.4.1. connection: Connect the metal hydroxide hydrogen reservoir to the gas instrument, and the pipe must be cleaned with hydrogen before the connection is completed

Way to prevent the air in the pipe from entering the system. The cleaning method of the pipeline is: open the metal hydroxide hydrogen reservoir valve and make the hydrogen flow through the pipeline into the atmosphere with moderate flow. In order to clean the pipeline, row it continuously for 5~10 seconds, and then tighten the connector and seal the system.

5.3.4.2. leakage detection: Check the interface with soapy water to ensure no air leakage.

5.3.4.3. hydrogen discharge: open the metal hydroxide hydrogen reservoir valve to allow for hydrogen discharge. In the process of hydrogenation of metal hydroxide hydrogen storage, the internal temperature is reduced, so the product needs to support a heating system, such as constant temperature and water bath, to ensure the smooth hydrogen discharge of metal hydroxide hydrogen storage, and can also heat up appropriately to improve the speed of hydrogen discharge. In general, 25℃ of water bath can delegate hydrogen can only reach 40% of the capacity, and 35℃ of water bath can reach 80% of the capacity, it is recommended to lower hydrogen in 35℃ of water bath. If the customer has higher requirements for hydrogen discharge speed and hydrogen discharge amount, the water bath temperature can be raised appropriately. The recommended water bath temperature is 50-60℃, and the hydrogen discharge amount can reach the nominal capacity.

#### **5.3.5 所 requires accessories**

Pressure relief valve for hydrogen charge: inlet pressure 15MPa, outlet pressure 6MPa;

Pressure relief valve for hydrogen discharge: inlet pressure of 3MPa, outlet pressure selected as required;

Connecting line: 3 Stainless steel pipe;

Heat exchanger device: water bath box.

Note: Regular testing and regeneration

When the product runs for a long time, leakage regularly, especially the connection and pipeline, to prevent hydrogen leakage; the capacity of the product (3 to 6 months), if the capacity is reduced to less than 80% of the rated capacity, regenerated.

Regeneration method: After the hydrogen reservoir is completely discharged, vacuum for a long time, and continue to vacuum under heating conditions, discharge impurities, water vapor, etc.: and then the hydrogen reservoir can be filled with hydrogen according to the hydrogen filling step.

#### **Transportation and storage conditions**

The metal hydrogen storage shall avoid extrusion during transportation and long-term storage.

## **6 Built-in lithium battery charging and discharge use methods and**

### **precautions**

#### **Note for 6.1.1 聚 compound lithium ion battery**

##### **6.1.1. charge**

Charging current and charging voltage shall not exceed the calibration value of the battery. If it exceeds the specified value, it may cause damage to the charging and discharge performance, mechanical performance and safety of the cell, and entry may lead to heat and leakage.

Battery charger must be capable of constant current and constant pressure charging, and the original charger is recommended;

The charging current of the single battery must be below 25.2V;

The temperature range is from 0 to + 45 °C when charging;

The voltage cannot exceed 26 V when charging.

##### **6.1.2. discharge**

The discharge current shall not exceed the following standards, and the discharge must be conducted within the scope of this standard.

The battery discharge current must be less than 15A;

The temperature range at discharge is from -20 to + 60 °C;

Battery discharge termination voltage shall not be less than 19.5V.

##### **6.1.3. over-discharge**

It should be noted that during the cell is not in use, it may use its self-discharge properties in some over-discharge state. To prevent overdischarge, the cell should be charged regularly to maintain its voltage above 23.0V. Excessive discharge can lead to the loss of cell performance and battery function.

#### **Polymer lithium-ion battery storage**

The long-term storage environment of the battery is: temperature -20~ + 35 °C

Relative humidity is 45 to 75%

During the battery storage period of nearly 1 year, the battery should be charged by 10% to 10% 50%.

#### **Polymer lithium-ion battery transport**

The battery shall be transported in a 10% to 50% charging state.

### **6.2 Polymer lithium-ion battery precautions:**

**6.2.1.** To prevent possible battery leakage, heat, or explosion, note the following precautions:

Do not removal the removed under any circumstances.

Do not immerse batteries in water or sea water.

Do not use or place batteries next to heat sources, such as fires, heaters, etc.

Do not heat or drop the battery into the fire.

Direct welding of batteries is prohibited.

No charging in a fire or hot environment.

Do not place the battery in a microwave or high pressure container.

Use or place batteries in high temperatures (e. g., strong sunlight or hot cars), otherwise it will cause overheating, fire, or function decline, and life to be reduced.

The 6.2.2. polymer lithium-ion battery does not theoretically have a flowing electrolyte, but in case of electrolyte leakage and contact with the skin, eyes or other parts of the body, it should immediately flush the electrolyte with clear water and seek medical treatment.

6.2.3. prohibit the use of damaged cells (cell plastic seal damage, damaged shell, smell of electrolyte gas, electrolyte leakage, etc.).

Batteries with electrolyte leakage or electrolyte odor should be kept away from the source of fire to avoid fire or explosion.

## 7 Installation method and precautions of heat companion sampling gun

### 7.1 Sampling gun parameters:

voltage:24V

power:110W

Length of heat pipe: 1500mm

Sampling steel pipe: 50mm


Temperature range: room temperature + 6- -200℃

### 7.2 Schematic diagram of the sampling gun

**1 Sampling end link handle 2 companion hot hose 3 inlet sample end link handle**

### 7.3 Schematic diagram of the installation link handle at the sampling end

- 1 Metal Filter 2 1 / 4 front sealing ring 3 1 / 4 rear seal gasket 4 M12\*1.25 stainless steel cap  
5 Inert stainless steel pipe 6 M12 \* 1 stainless steel cap 7 1 / 4 rear sealing gasket 8 1 / 4 front sealing ring  
9 Steel pipe link head 10 seals the O-ring 11 microhole filter film 12 mesh plate 13 sampling end handle

### 7.4 Schematic diagram of the links to the sampling gun

The sampling gun into the sample end of the joint red dot aligned with the sampling gun interface base red card slot inserted vertically when the sound sound link is in place.

### 7.5 Pull of sampling gun

Hold the upper part of the red point of the link head of the sampling gun and pull up vertically to pull out.

## 8 Software operation instructions

### install

This section describes the installation of the software. The legitimate HSNMHC400 special analysis software, such as operating system, software, data damage needs to be reinstalled, please follow the following steps to install.

#### 8.1 8 method of erection

This section describes the installation of the software.

There are 2 steps to use the software from installation to use:

##### 8.1.1 install

The 8.1.2 序 column number registers

### 8.2 Before installation

Before installing the software, check the following:

The operating system running this software is Windows 10. Please confirm the operating system of the PC.

Prepare the registration number serial code (in the U disk or instructions).

#### 8.2.1

Insert the software U disk to find the (HSNMHC400 software installation package) compressed file

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This software is green software without installing decompression, after decompression can find the HSNMHC400.exe executable to create a desktop shortcut.

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### 8.2.2 软件 Part registration

The first operation of the software requires registration code registration before operation

Double-click on the HSNMHC400 software shortcut icon

The first operation needs to input the registration code: in this manual and software U disk, or ask the installation after-sales engineer.

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Enter the registration code to determine it.

## 8.2 Software operation method

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Manual editing of the computer IP address is segment 192.168.1.100

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IPv4

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192.168.1.100

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1

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
IPv6

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### 8.3 Software interface introduction:

| number | explain                                                                                                    |
|--------|------------------------------------------------------------------------------------------------------------|
| 1      | (Function toolbar) Display: analysis, setting, calibration, spectral map query, historical data, log query |
| 2      | Real-time analysis of the data display regions                                                             |
| 3      | Run the control area                                                                                       |
| 4      | Temperature setting control zone                                                                           |
| 5      | Real-time signal display area                                                                              |
| 6      | Operating status indication area                                                                           |
| 7      | Battery power display area                                                                                 |
| 8      | Exit software                                                                                              |

9 Software operation function introduction

9.1 analyse

9.1.1 click the **关闭吹扫** button to sweep the sampling gun and residual pollutants in the sampling pipe  
**打开吹扫** to turn off the purpower  
Can.

9.1.2 click [Start Analysis], the device performs the analysis program, stop by itself after the program, or click [Stop Analysis] to terminate the analysis program in advance


The 9.1.3 temperature setting

Enter the setting temperature to check the heating part rear point [set] key to download the parameters to the equipment point [one key start] to perform the heating operation.

#### 9.1.4 figure adjustment

Press the left mouse button in the map display area, move the mouse can move the map position, the left key circle to enlarge the area to zoom in, click the mouse right click [restore the default size] to restore the original state

#### The 9.1.5 status indicator bar

The initial state is blue and green is ready, reaching the instrument setting parameter value.

[Total number] is expressed as the total number of analysis for this analysis


[Current times] represents the current number of execution

[Single-process timing] represents the total analysis time of a single cycle

[Time] Display the current time


9.1.6 点 ignition automatically 10 minutes after default temperature control green. If manual ignition, click [Ignition] button [OK]

9.2 set up


Point [Settings] enter the initial password 123456 [confirm] to enter the setting interface

9.2.1degree setting

Temperature calibration parameters, follow the steps:

Prepare the calibrated mercury thermometer, insert the parts to calibrate the temperature, read the mercury thermometer temperature value after the data stabilizes, fill in the [real temperature] column, set the temperature, fill in the [set temperature] column, and click [Download the temperature calibration parameters].

### 9.2.2 detector

[Ignition] is equivalent to the ignition function in the analysis interface, and the manual ignition switch.

[Fire threshold] This parameter is to determine the ignition success value of the device. The ignition is considered successful when the baseline is above the set point. If lower than the setting value and judged as fire extinguishing.

[Ignition status] Fire extinguishing is OFF fire NO

[Ignition times] When the equipment determines that the fire extinguishing state is the number of automatic ignition performed, the maximum is 5 times

[Amplification] Sensor sensitivity adjustment

[Phase] Adjust the peak direction

[Settings] All parameters change and need to click the Settings key to confirm

### 9.2.3part settings

[Sampling] The start time of the mining sample is generally 00:01

[Sampling] End time of off sampling, 00:20

[Analysis] The start time of the open analysis, which generally starts at 00:21

[Analysis] End time of turn-off analysis, time after methane production.

Set the breakpoint [Download the event parameter]

Cycle parameter settings

[Number of process cycles] Number of samples analyzed for an analysis process.

[Single analysis time] The total analysis time for a single process includes the sampling and analysis time.

[Single delay time] Generally set to the same time of the sampling pass

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192.168.0.105  
192.168.0.106  
192.168.0.100

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[Map file name setting] These three items are used for the control selection of the file after the end of the experiment.

The first item indicates that the file is kept at the date and time of the experiment and is saved automatically;

The second way is to save the file name you entered as the date and time of the experiment, and use automatic saving;

The third way is to save the file that you typed as the previous file name, with 000n ++ count, and automatically add one for each saved file system count

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[Settings] Unit conversion, ppm and mg/m<sup>3</sup> units can automatically convert the settings and click the [Save] button

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PPM

mg/m<sup>3</sup>

[Power calibration] Built-in lithium battery power display calibration setting, click [Save] button

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[Communication IP settings] Synchronize the local computer IP address segment, wifi link settings

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
| parameter                          | scope          | function declaration                                                                                                                                                                                                       |
|------------------------------------|----------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Minimum half peak width            | 0-50           | Removing peaks with small or large peak width ratio set values is like the screen based on the half height width of the minimum effective peak. The minimum half peak width can be estimated from the spectrograph display |
| Maximum half peak                  | 1000 more than |                                                                                                                                                                                                                            |
| gradient                           | 0-99999.9      | Peak detection sensitivity, used to determine the beginning and end point of the peak                                                                                                                                      |
| minimal area                       | 0-99999.9      | Peaks with the area (peak height) smaller than the set point were used to remove area peak height relative to smaller unrelated peaks (e. g., baseline noise) during the analysis                                          |
| The minimum peak is high           | 0-99999.9      |                                                                                                                                                                                                                            |
| The number of smooth               | 0-30           |                                                                                                                                                                                                                            |
| Remove the maximum (minimum) value | 0-99.9         | Remove a minimum content or maximum content value                                                                                                                                                                          |
| Number of window sampling points   | 0-30           |                                                                                                                                                                                                                            |

The 9.3.2 校准 quasi-curves are manually selected automatically

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
[Manual selection]: Four optional high concentration in the low concentration of a single point.

Select the response curve according to the approximate concentration point of the sample reagent and click [Save the calibration parameters]

[Automatic selection]: After point selection, set low and medium, high concentration peak area boundary parameters, such as medium and low concentration setting 100, medium and high concentration setting 200 indicate that the low concentration calibration 100-200 when the total hydrocarbon area is less than 100, and the high concentration calibration is higher than 200. Production of the

### 9.3.3quasi-curves


#### 9.3.3.1 selection [spectral calculation]

#### 9.3.3.2 [Open] spectral map file and select the corresponding concentration standard spectral

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
9.3.3.3 Click [Start peak search] to adjust the number of sampling points in the integral parameters, and the smoothing number parameters until the total hydrocarbon and methane peak time and area appear.


9.3.3.4 calibration in the low concentration curve, the medium to high concentration curve production method is not repeated.

Select [low concentration curve] to query the time of total hydrocarbon peak time and area in [spectral chart calculation] in [serial number] 1 [peak recognition time] column, [component name] enter name [concentration 1] enter the corresponding standard gas concentration value, and [area 1] enter peak area. Enter the time of methane in the column [No.] 2 [peak recognition

time], [component name] enter the name [concentration 1] enter the corresponding standard gas concentration value, and [area 1] input the peak area. So on, repeat 1-3 steps to retrieve different gas concentration files for calibration for a maximum of 10 points. The gas concentration document shall be selected following the principle from low to high.

9.4 [Map query] Select start and end time Click [Query] to retrieve map data for this time period.

9.5 [Historical data] Instrument model, instrument number, measurement place, measurement time. Check the point according to the actual situation [Save data]

9.5.1[creation time] Select the start and end time of the creation time. Click [Query] to count the average value of the total hydrocarbon, methane and non-methane total hydrocarbon concentrations during this time period. If a data is wrongly selected, click [Delete data] to delete it.


After connecting the intelligent smoke gun, Bluetooth can check the working condition data to be printed, such as: ambient temperature, pressure, etc

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9.6 [Log query] is divided into operation log alarm log operation log, select the start and end time of log time, and click [query] to query the equipment operation of this time period.Click [Delete] to delete this log content.

## 10 Analyze the fault judgment and treatment methods

| fault phenomenon                                             | Troubleshooting                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              | Inspection method and repair                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
|--------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>There is no spectral peak</b>                             | 12. The amplifier power supply is disconnected.<br>13. The ion line is broken.<br>14. No air no not over.<br>15. Poor contact with the integrator / chromatographic workstation signal line.<br>16. Integrator / chromatograph workstation failure<br>17. The sample inlet temperature was too low so that the sample was not vaporized.<br>18. Microinjectors were blocked.<br>19. Air leakage of silicone pad of sample inlet.<br>20. The chromatcolumn installation connection is not standard.<br>21. The fire is not lit (FID).<br>22. The FID polarization voltage is unconnected or poorly contacted. | 12. Check the amplifier, the fuse.<br>13. Check the ion line.<br>14. Check that the flow path is blocked and that the gas supply is exhausted.<br>15. Check the integrator / chromatographic workstation wiring.<br>16. Trouthe integrator / chromatograph workstation.<br>17. Increase the sample inlet temperature.<br>18. Replace the syringe.<br>19. Replace the sampler silicone pad.<br>20. Reinstall and tighten the chromatcolumn.<br>21. re-light.<br>22. Connect the polarized voltage to exclude poor contact. |
| <b>The tention time was normal but sensitivity decreased</b> | 8. The decay is too large<br>9. Not a sufficient sample size<br>10. Loss during sample entry<br>11. Injector leakage or blockage<br>12. Air leakage<br>13. Improper selection of Hydrogen and Air Flow (FID)<br>14. Unnormal polarization voltage (FID)                                                                                                                                                                                                                                                                                                                                                      | 8. Reduce the attenuation and increase the scale range<br>9. Increase sample intake<br>10. Ensure that all the samples go into the system<br>11. Replace the syringe<br>12. Detect the air leakage and eliminate the air leakage<br>13. Adjust the hydrogen gas and the air flow rate<br>14. Check and troubleshoot the polarization voltage fault                                                                                                                                                                        |
| <b>tailing peak</b>                                          | 5. The intake temperature is too low<br>6. Vaporization pipe pollution of sample inlet or residual sample inlet pad                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          | 6. Readjust the sample inlet temperature<br>7. Clean the inlet or remove the                                                                                                                                                                                                                                                                                                                                                                                                                                              |

|                                                         |                                                                                                                                                                                                  |                                                                                                                                                                                             |
|---------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                                                         | <p>7. The chromatcolumn temperature is too low</p> <p>8. The feeding technique is too poor</p> <p>The 5. chromatcolumn was improperly selected</p>                                               | <p>residue from the inlet pad</p> <p>8. Increase the chromatcolumn temperature</p> <p>9. Improve the sample entry technology</p> <p>10. Select the appropriate chromatcolumn</p>            |
| <b>Extend the tongue peak</b>                           | <p>3. The sample size was too large for the chromatcolumn overload</p> <p>4. Samples are condensed in the system</p>                                                                             | <p>3. Reduce sample size</p> <p>4. Increase the temperature and aging system</p>                                                                                                            |
| <b>Chmatographic peaks were not separated</b>           | <p>4. The chromatcolumn temperature is too high</p> <p>5. The chromatography column is too short</p> <p>6. The gas flow rate is too high</p> <p>The 4. chromatcolumn was improperly selected</p> | <p>5. Reduce column temperature</p> <p>6. Select longer chromatcolumns</p> <p>7. Reduce the gas flow rate</p> <p>8. Select the appropriate chromatcolumn</p>                                |
| <b>flat peak</b>                                        | <p>3. The amplifier is saturated</p> <p>4. Integrator / chromatograph workstation failure</p>                                                                                                    | <p>3. Reduce the sample size or reduce the amplifier sensitivity</p> <p>4. Trouthe integrator / chromatograph workstation</p>                                                               |
| <b>There are irregular fluctuations at the baseline</b> | <p>5. The instrument is not down to earth</p> <p>6. Detector contamination</p> <p>7. The airflow selection is improper</p> <p>8. The amplifier is not good</p>                                   | <p>5. Improve grounding</p> <p>6. Clean the detector</p> <p>7. Select the appropriate gas flow rate</p> <p>8. Check the amplifier</p>                                                       |
| <b>The baseline noise was large</b>                     | <p>6. Chromatic column pollution</p> <p>7. The gas is not pure</p> <p>8. bad earth</p> <p>9. Detector contamination</p> <p>10. Poor detector cable contact</p>                                   | <p>6. Replace the color column</p> <p>7. Purification of the gas meets the requirements</p> <p>8. Improve grounding</p> <p>9. Clean the detector</p> <p>10. Replace or repair the cable</p> |
| <b>Baseline single-direct ional drift</b>               | <p>3. The gas is gradually over</p> <p>4. System air leakage</p>                                                                                                                                 | <p>3. Replace the cylinder</p> <p>4. leak detection</p>                                                                                                                                     |

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